FITPADDY®

FITMASK

Maschera LED viso / LED Facial mask
Manuale d'uso / User Manual

IT – Italiano
EN – English
FR – Français
ES – Español
PT – Português
DE – Deutsch
NL – Nederlands
PL – Polski
CS – Čeština
SK – Slovenčina
HU – Magyar
BG – Български
GR – Ελληνικά
SR – Srpski
HR – Hrvatski
SL – Slovenščina
SE – Svenska
FI – Suomi
DK – Dansk
NO – Norsk
RO – Română
LT – Lietuvių
LV – Latviešu
ET – Eesti
UK – Українська
IS – Íslenska
TR – Türkçe
RU – Русский
AL – Shqip

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INTENDED USE
FITMASK is an LED facial mask designed for your cosmetic skincare routine. Its four light modes help visibly improve skin appearance. The product is intended exclusively for indoor domestic use and is not a medical device.

PACKAGE CONTENTS
• 1 FITMASK • 1 Controller • 1 USB-C connection cable • 1 USB charge cable
• 1 Adjustable headband • 2 Eye shields • 1 User manual

PREPARATION AND USAGE
• Fully charge the controller before first use using a 5V⎓ 1A USB power adapter.
• Check the product and cables for any signs of damage.
• Clean and dry your face, remove glasses, contact lenses, and false eyelashes, then put on the eye shields.
• Connect the mask to the controller, put it on, and adjust the headband.
• Press and hold the ON/OFF button to turn the device on.
• Briefly press ON/OFF to cycle through the 4 light modes, and press the M button to set the timer (from 10 to 30 minutes).
• Keep your eyes closed while using the device.
• The device turns off automatically when the set time expires.
• We recommend 3 to 5 sessions per week.

MODES
• Red Mode (Red + IR): Combines red and infrared light to help make skin look visibly firmer and more radiant.
• Blue Mode (Blue + IR): Combines blue and infrared light to help improve the appearance of blemish-prone skin.
• Yellow Mode (Red + Yellow + IR): Combines red, yellow, and infrared light to help give the skin a more even and glowing look. During use, the visible light may appear predominantly red.
• Purple Mode (Red + Blue + IR): Combines red, blue, and infrared light for a complete cosmetic routine.

CHARGING
Only use the supplied cables and a 5V⎓ 1A USB power adapter. Do not use the device while it is charging. Unplug the controller from the power source once charging is complete.

WARNINGS
• Use exclusively for the intended purpose.
• Keep your eyes closed while using the device and use the provided eye shields.
• Do not look directly at the active LEDs.
• Do not use if you are photosensitive, undergoing photosensitizing medication therapy, or following invasive aesthetic treatments without consulting a doctor.
• Consult your doctor if you are pregnant, have epilepsy, a pacemaker, implantable devices, or other medical conditions.
• Stop using the device immediately if you experience pain, irritation, or any other adverse effects.
• Do not use while charging, with wet hands, or if the device comes into contact with water.
• Do not submerge, dismantle, modify, or use the product if damaged.
• Do not expose to fire, heat sources, or high temperatures.
• Keep out of reach of children.

CARE AND MAINTENANCE
Clean the mask with a soft, slightly damp cloth. Do not use solvents or harsh detergents. Store in a cool, dry place.

DISPOSAL
Do not dispose of the product with household waste. Take it to an authorized WEEE (RAEE) collection center.

SPECIFICATIONS
• Model: FITMASK
• Power Supply: 5V DC 1A
• Battery: Li-ion 3,7 V – 2100 mAh (7,77 Wh)
• Rated Power: 8W
• LEDs: 96
• Charging Port: USB-C
• Timer: 10–30 minutes
• Modes: 4
• Autonomie : environ 80 minutes
• Usage domestique en intérieur

SYMBOLS AND COMPLIANCE
• CE: The product complies with the essential requirements of applicable European Union legislation.
• WEEE (RAEE): The product must not be disposed of with regular household waste, but must be taken to an authorized collection center for the separate recycling of electrical and electronic equipment.
• EU Declaration of Conformity: FITPADDY S.r.l. declares that the FITMASK product complies with the requirements of Directives 2014/30/EU (EMC) and 2011/65/EU (RoHS). The full EU Declaration of Conformity is available at: https://fitpaddy.com/declaration-of-conformity/FITMASK

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
FITPADDY®

User Manual
FITMASK
LED Face Mask

Product overview
FITMASK is an LED face mask designed for cosmetic skincare and facial wellness. It uses different light wavelengths to provide a non-invasive cosmetic treatment, without chemicals, suitable for home use. The different light modes allow treatment to be customised according to individual cosmetic needs, helping improve the skin's overall appearance and support a daily skincare routine. The device is intended solely for cosmetic purposes and is not a medical device. It is not intended to diagnose, treat, cure or prevent disease.

Product features

Effective
Equipped with 96 LEDs, the mask uses different light wavelengths for cosmetic treatment of facial skin.

Light modes
The 4 light modes allow the cosmetic treatment to be customised to suit individual needs and preferences.

Safety
Non-invasive cosmetic treatment designed for home use in accordance with the instructions for use.

Convenient
Covers the entire face in a single application and features an automatic timer for simple, practical use.

Comfort
Made of soft silicone, it fits comfortably over the face and helps reduce external light entering the mask.

Customisable
The different treatment modes allow the routine to be customised according to individual preferences.

Practical
Controller with intuitive controls for easy management of the modes and timer.

How to use
It is recommended to use the device 3 to 5 times per week.

Clean the skin thoroughly and gently pat it dry with a soft towel before using the device.

Insert the elastic strap into the mask.

Connect the mask's USB Type-C connector to the controller's USB Type-C port.

Open the mask, place it on the face and secure it around the back of the head using the elastic strap. Adjust the strap for a comfortable fit and make sure the eye protectors are correctly positioned. Do not look directly at the light emitted by the LEDs.

Press and hold the power button for approximately 1.5 seconds to switch on the device. The default mode at start-up is red light. Briefly press the power button to select the desired light mode.

Briefly press the M button to select the treatment duration. Five settings are available: 10, 15, 20, 25 and 30 minutes.

At the end of the selected treatment time, the device switches off automatically. To switch it off manually at any time, press and hold the power button for approximately 1.5 seconds.

1. Anti-ageing Mode (R + NIR)

Light emitted:
- Red light (peak wavelength 633 nm)
- Near-infrared light (NIR, peak wavelength 830 nm)

Description:
This mode combines red light and near-infrared light (NIR) for a cosmetic skin treatment, helping improve its appearance and promote a more even, radiant-looking complexion.

2. Target Mode (B + NIR)

Light emitted:
- Blue light (peak wavelength 415 nm)
- Near-infrared light (NIR, peak wavelength 830 nm)

Description:
This mode combines blue light and near-infrared light (NIR) and is intended for the cosmetic treatment of blemish-prone skin, helping improve its appearance.

3. Recovery Mode (R + Y + NIR)

Light emitted:
- Red light (peak wavelength 633 nm)
- Yellow light (peak wavelength 590 nm)
- Near-infrared light (NIR, peak wavelength 830 nm)

Description:
This mode combines red, yellow and near-infrared light (NIR) for a gentle cosmetic treatment, helping improve the skin's appearance and promote a brighter, more even-looking complexion.

4. Balance Mode (R + B + NIR)

Light emitted:
- Red light (peak wavelength 633 nm)
- Blue light (peak wavelength 415 nm)
- Near-infrared light (NIR, peak wavelength 830 nm)

Description:
This mode combines red, blue and near-infrared light (NIR) for a complete cosmetic facial treatment. It is designed to provide a balanced treatment experience and help improve the skin's overall appearance.