FITPADDY

MANUALE D'ISTRUZIONI / USER MANUAL

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INTENDED USE
FITROLLER is a roller muscle massager designed for massaging the abdomen, back, and legs. It is intended for home use and to promote muscle relaxation after physical activity or during recovery periods.

SAFETY GUIDELINES
- Use exclusively for muscle massage and relaxation.
- Keep out of reach of children.
- Do not expose to extreme temperatures, humidity, or direct sunlight.
- Avoid impacts, drops, or excessive pressure on the device.
- Do not disassemble or modify the product.
- Stop use immediately in case of pain, discomfort, or abnormal reactions.
- Do not use if you have hypertension, heart disease, thrombosis, tumors, osteoporosis, or acute illnesses.
- Not suitable for use by children or people with reduced mobility.

OPERATION
- Press and hold the power button for 2 seconds to turn on FITROLLER.
- When switched on, the device automatically starts at level 1.
- A single press of the button increases the speed (levels 1 to 6).
- Quickly press the button twice to change the rotation direction.
- Press and hold the button for 2 seconds to turn off the device.
- FITROLLER features automatic shut-off after 30 minutes of continuous use.
- Use on dry skin without applying excessive pressure.

CHARGING
- Charge FITROLLER using a USB Type-C cable.
- Charging time: approximately 5 hours.
- Do not use the device while charging.
- Use certified chargers only.

WARNINGS
- Do not use if the device shows visible damage.
- Do not cover the device during operation to prevent overheating.
- Do not use outdoors.
- After use, allow the device to cool down before storing it.

CARE AND MAINTENANCE
- Clean with a dry or slightly damp cloth.
- Do not use aggressive detergents, solvents, or abrasive products.
- Do not immerse in water.
- Store in a cool, dry place away from heat sources.

LIMITATION OF LIABILITY
FITPADDY shall not be held responsible for damage or malfunctions resulting from improper use, unauthorized modifications, impacts, drops, exposure to extreme conditions, a discharged battery, or inadequate maintenance.

DISPOSAL
Dispose of the product in accordance with local regulations on electronic waste (WEEE). Do not dispose of with household waste.

TECHNICAL SPECIFICATIONS
- Manufacturer code: CM2529
- Power supply: Rechargeable lithium battery
- Voltage: 7.4V
- Battery capacity: 4000 mAh
- Speed levels: 6
- Battery life: up to 15 hours (level 1, no load)
- Charging time: approximately 5 hours
- Automatic shut-off: 30 minutes
- Dimensions: 525 × 140 × 70 mm
- Net weight: 1341 g

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SERVICE CLIENT
Site Web: www.fitpaddy.com
Email: customercare@fitpaddy.com
Fabricant: FITPADDY Srl
Via dell'Artigianato 16, 20835 Muggiò (MB), Italie

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EN - This product complies with Directive RoHS 2011/65/EU and 2015/863/EU.
The product bears the CE mark and complies with the safety standards set by the European Union.
FITPADDY Srl declares that the product complies with Directive 2014/53/EU.
The full text of the EU Declaration of Conformity is available at the following Internet address: https://fitpaddy.com/declaration-of-conformity/FITROLLER

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FITPADDY Srl shpall se produkti është në përputhje me direktivën 2014/53/UE.
FITPADDY Srl deklaruoja, kad produktas atitinka 2014/53/UE direktyvą.
FITPADDY Srl paziņo, ka produkts atbilst direktīvai 2014/53/UE.
FITPADDY Srl kuulutab, et toode vastab direktiivile 2014/53/UE.
FITPADDY Srl prohlašuje, že výrobek je v souladu s nařízením 2014/53/EU.
FITPADDY Srl vyhlasuje, že výrobok je v súlade so smernicou 2014/53/EÚ.